

Introduction to G-Flash™, an Alternative to Traditional Flash Chromatography



Kenneth J. James, Ph.D.


G-Flash Chromatography, LLC

Abstract

Traditional peptide purification relies on solvent-intensive Flash Chromatography, posing significant environmental and safety challenges. This study evaluates G-Flash™, a supercritical CO₂ (sc-CO₂) platform, as a green alternative for purifying the synthetic peptide VGVAPG-COOH. By comparing N-terminal Fmoc-protection against the use of chaotropic mobile phase modifiers, we demonstrate that masking peptide polarity improves peak resolution.

Inquiry Justification


First, current peptide purification relies on hazardous organic solvents, necessitating a transition toward sc-CO₂ to improve laboratory safety and reduce environmental costs. Second, while traditional Supercritical Fluid Chromatography is often too complex and expensive for routine use, the G-Flash™ system provides a streamlined, cost effective alternative for green chromatography. Finally, this study addresses the technical challenge of separating polar peptides in non-polar mobile phases by evaluating whether N-terminal Fmoc protection or chaotropic modifiers offer a more effective path to purification.




Research Question

Can the purification of synthesized peptides be achieved more effectively in a supercritical CO₂ (G-Flash™) environment by using N-terminal Fmoc-protection to reduce polarity, compared to the use of chaotropic mobile phase modifiers?

Background and Data



Unprotected Peptide: 220 nm, 80-60% scCO₂ gradient, 45 ml/min Total Flow; C18 Spherical Silica Cartridge



Fmoc Protected Peptide: 220 nm, 90-60% scCO₂ gradient, 50 mL/min Total Flow; C18 Spherical Silica Cartridge

Results

Peptide purification was evaluated in a scCO₂ environment, comparing N-terminal Fmoc-protection versus Hofmeister salt modifiers. Optimization of gradients and cosolvents achieved baseline separation for both compounds, verified by blank injections of DMSO, TFA, and NH₄OH. The unprotected peptide reached nearly 70% peak purity with a resolution factor, R_s of 2.0, with 3 additional peaks as synthesis residual. In contrast, the Fmoc-protected peptide demonstrated superior resolution, $R_s=3.1$, and improved isolation from solvent peaks with some residual impurities in the DMSO peak. These results indicate that reducing analyte polarity through Fmoc-protection is more effective than traditional modifiers for maximizing chromatographic clarity in scCO₂ systems.



Fig. 1 Sorbtech Cartridges

Top Cartridge (Blue):
Normal Phase 70Å, 20-45 um, 12g
Bottom Cartridge (Purple):
C18 Reverse Phase 100Å,
20-45 um, 20g

Experimental Design:
2 Factor Runs x 2 levels x 2
Replicates = 8 Total Runs

Run	Cosolvent	Cartridge	Gradient
1	Cosolvent A	Normal Phase	Gradient A
2	Cosolvent B	Normal Phase	Gradient A
3	Cosolvent A	Reverse Phase	Gradient A
4	Cosolvent B	Reverse Phase	Gradient A
5	Cosolvent A	Normal Phase	Gradient B
6	Cosolvent B	Normal Phase	Gradient B
7	Cosolvent A	Reverse Phase	Gradient B
8	Cosolvent B	Reverse Phase	Gradient B

Fig. 2: Outline of Factorial Design Experiment

Conclusion

This study demonstrates that N-terminal Fmoc-protection has strong potential for peptide purification in supercritical CO₂ (G-Flash™) environments compared to the use of chaotropic modifiers. While Hofmeister salts facilitated separation for the unprotected peptide, $R_s = 2.0$, the resulting chromatogram naturally reflected the 70% purity of the starting material. In contrast, the Fmoc-protected peptide achieved a higher resolution factor of 3.1, more effectively isolating the target analyte from both the DMSO solvent and inherent compound impurities following the peak at 1 minute. These results indicate that reducing analyte polarity through structural protection provides a more robust and high-resolution chromatographic profile, allowing for more effective purification in scCO₂ systems than ionic mobile phase additives.

Meet the Author

Kenneth J. James, PhD



Ken earned his Ph.D from the University of Delaware, where he specialized in analytical spectroscopy and co-developed the commercialized “Brill Cell” for material decomposition. A seasoned entrepreneur, he founded

Supercritical Fluid Technologies and G-Flash, LLC, while authoring several key patents in chromatography. He now leads all technical development for G-Flash, LLC, focusing on the commercialization of Green Flash Chromatography (GFC).

References

[1] Regalado E.L., et al. “Chaotropic Effects in Sub/Supercritical Fluid Chromatography via Ammonium Hydroxide in Water-Rich Modifiers: Enabling Separation of Peptides and Highly Polar Pharmaceuticals at the Preparative Scale.” *Analytical Chemistry*, vol. 91, no. 21, 2019, pp. 13907–13915., <https://doi.org/10.1021/acs.analchem.9b03408>.

Contact Us

Kenneth J. James Ph.D.
President and Director of Technology
550 S. College Ave
Newark, DE 19713



+1 (302)-738-3420 x201
Ken.James@GreenFlashChromatography.com